

Infant Incubators



NEO Shield- Series

Creating New Dimensions for Neonatal & Child Care



- Hi-end NICU Incubator
- Micro Computer Servo Controlled
- Baby skin mode & Air mode
- Hi-Temperature cut off facility
- Optional online baby weighing scale

Microprocessor based Servo controlled Temperature system.

Control Modes: Air mode and Baby Skin Mode.

Humidity control system with humidity concentration display system.

Set Temp, Air Temp, Baby Skin Temp, Timer and Heating power are displayed separately.

Self-testing function, various failure alarms by audible and visual.

Embedded integrated sensor module.

Triple protection for over temperature with separate cut off device for extra safety system.

Double wall hood, 4 operating windows and 2 iris ports with Independent Locking device for front door



Manufactured by:



Aster Medical Systems

Phone : +91 9030678639

Email : astermedicalsyste.ms@gmail.com

Website : www.astermedicalsyste.ms.com



NEOKRAFT™
ISO 13485 Certified Company