

Infant Resuscitator



RESIPUFF 1050

Infant T-Piece resuscitator
Lightweight easy to operate
No power or battery is required



Resipuff 1050 Infant T-Piece resuscitator is a lightweight stand alone unit suitable for delivery wards, Neonatal care units, Operating rooms and Infant transport resuscitation purpose.

- Infant resuscitator will safely inflate the baby's lungs and provide optimum oxygenation by delivering consistent PIP & PEEP with each breath avoiding the risks associated with under or over inflation of uncontrolled pressures
- The patient T-Piece can connect to Neonatal masks or Endotracheal tubes directly.
- Different sizes of face mask which provides flexible in usage for different babies
- Inspiratory high pressure safety cut-off relief setting
- Peak inspiratory pressure setting (PIP)
- Positive end expiratory pressure setting (PEEP)
- Works on low flow from 5 L/m to 15 L/m
- Medical grade monometer to monitor PIP & PEEP delivered



NEOKRAFTTM
ISO 13485 Certified Company